

Closure Accessory and Supply Ordering information

Supply Item	Description/Comments	Ordering Information
VNUS® Closure® Procedure Packs	Pre-packaged kits with Closure procedure accessories	
Pack with no introducer sheath set	<ul style="list-style-type: none"> • Sterile • Peel pouch 	Part #: CPP VNUS Medical Technologies 888-797-8346
Pack with 6F introducer sheath set	<ul style="list-style-type: none"> • 6F 7cm introducer sheath set with .035" 45cm DE guidewire • Sterile • Peel pouch 	Part #: CPP-6F VNUS Medical Technologies 888-797-8346
Pack with 8F introducer sheath set	<ul style="list-style-type: none"> • 8F 7cm introducer sheath set with .035" 45cm DE guidewire • Sterile • Peel pouch 	Part #: CPP-8F VNUS Medical Technologies 888-797-8346
Introducer Sheaths Sets		
VNUS 6F Introducer Sheath Set, 7cm length	<ul style="list-style-type: none"> • 6F 7cm introducer sheath set w/ .035" 45cm DE guidewire 	Part #: IS-6F07 VNUS Medical Technologies 888-797-8346
VNUS 8F Introducer Sheath Set, 7cm length	<ul style="list-style-type: none"> • 8F 7cm introducer sheath set w/ .035" 45cm DE guidewire 	Part #: IS-8F08 VNUS Medical Technologies 888-797-8346
VNUS 6F Introducer Sheath Set, 11cm length	<ul style="list-style-type: none"> • 6F 11cm introducer sheath set w/ .035" 45cm DE guidewire 	Part #: IS-6F11 VNUS Medical Technologies 888-797-8346
VNUS 8F Introducer Sheath Set, 11cm length	<ul style="list-style-type: none"> • 8F 11cm introducer sheath set w/ .035" 45cm DE guidewire 	Part #: IS-8F11 VNUS Medical Technologies 888-797-8346
Daig Fast-Cath	6F or 8F introducer set <ul style="list-style-type: none"> • Set includes 5cm sheath, guidewire, and dilator • Requires 18g (or larger), thin-wall needle to accept a 0.038" guide wire. 	Part #: 6F: 406164 8F: 406170 Daig Corp, St. Jude Medical Minnetonka, MN 800-328-3873
J & J Cordis Avanti	6F or 8F introducer set <ul style="list-style-type: none"> • Consists of 11cm sheath, guidewire, and dilator • Requires an 18 gauge (or larger), thin-wall needle to accept an 0.038" guide wire. 	Part #: 6F: 504-606X 8F: 504-608X Cordis Corp. Miami, Florida 800-327-7714
Micropuncture Kits		
VNUS 4F Micropuncture Kit	<ul style="list-style-type: none"> • 4F, 7cm micropuncture kit w/ dilator, .018" 45cm guidewire, 21g access needle • Use to size up to a 6 or 8F sheath 	Part #: MPK-4F VNUS Medical Technologies 888-797-8346
VNUS 5F Micropuncture Kit	<ul style="list-style-type: none"> • 5F, 7cm micropuncture kit w/ dilator, .018" 45cm guidewire, 21g access needle • Use to size up to a 6F or 8F sheath 	Part #: MPK-5F VNUS Medical Technologies 888-797-8346
Supply Item	Description/Comments	Ordering Information

Micropuncture Kits		
Cook	<ul style="list-style-type: none"> • 4F, 10cm micropuncture kit w/ dilator, .018" 45cm guidewire, 21g access needle • Use to size up to a 6F or 8F sheath 	Part #: 216347 Cook Corp. Bloomington, Illinois 800-457-4500
Micro Introducer Set		
VNUS 6F Micro Introducer Set	<ul style="list-style-type: none"> • Combines the features of a micropuncture kit and standard introducer sheath • 6F, 7cm sheath w/ dilator tapering to 2F, .018" 45cm guidewire, 21g access needle 	Part #: MIS-6F07 VNUS Medical Technologies 888-797-8346
VNUS PTFE-Coated Guidewires	<ul style="list-style-type: none"> • .025", 150cm length • PTFE (Teflon)-coated • Double-ended: 1.5mm J and straight • Used to help navigate the Closure <i>PLUS</i> catheter 	Part #: GW-025-150-DE-1.5J VNUS Medical Technologies 888-797-8346
VNUS PTFE-Coated Guidewires	<ul style="list-style-type: none"> • .025", 260cm length • PTFE (Teflon)-coated • Double-ended: 1.5mm J and straight • Used to help navigate the Closure <i>PLUS</i> catheter 	Part#: GW260 VNUS Medical Technologies 888-797-8346
Boston Scientific/Microvasive	<ul style="list-style-type: none"> • PTFE (Teflon)-coated • 150cm length • .025" 	Part #: 620-107 (regular shaft) or 620-108 (stiff-bodied) Boston Scientific N. Quincy, MA 800-225-3226
Boston Scientific/ Microvasive	<ul style="list-style-type: none"> • Glidewire • 150 and 260cm lengths • .025" 	Part #: 630-093 (150 cm) or 5199 (260 cm) 800-225-3226
Y-Connector		
B. Braun Medical, Inc.	<ul style="list-style-type: none"> • 9.5F lumen • Single Y-connector w/ plastic guidewire insertion tool 	Part #: 610 402 B. Braun Medical, Inc Bethlehem, PA 800-854-6851
Esmark Wraps		
Bioseal	<ul style="list-style-type: none"> • Latex-free • 6" by 4 yards • Sterile 	Part #: 04964/20 Bioseal Placentia, CA 714-528-4695
Esmark Wraps	•	
DeRoyal Industries	<ul style="list-style-type: none"> • Latex-free • 6" by 4 yards • Sterile • CE mark • 20 per case 	Part #: 35-612 DeRoyal Industries Lake Powell, TN 800-251-9864
Supply Item	Description/Comments	Ordering Information
Needles		

VNUS Percutaneous Entry Needle	<ul style="list-style-type: none"> • 19g UTW (ultra-thin wall) • Accommodates up to a .035" guidewire 	Part #: AN-19G VNUS Medical Technologies 888-797-8346
Cook	<ul style="list-style-type: none"> • 19g UTW • Accommodates up to a 0.035" guidewire 	Part #: BSDN-19UT-7.0 Cook Corp. Bloomington, Illinois 800-457-4500
Echo-Coat Ultrasound Needles	<ul style="list-style-type: none"> • 22 gauge 3.5" spinal needle • Echolucent coating allows easy visualization under duplex ultrasound • For ultrasound-guided tumescent infiltration 	Part #: EC-SP22G09CM STS Biopolymers, Inc. Henrietta, NY 716-321-1130
Transducer Cover		
GE Sterile Intraoperative Probe and Cord Cover, with or without gel	<ul style="list-style-type: none"> • 6" tapered to 3" width x 96" length • Latex-free polyethylene • CE marked • 12 covers/box 	Part #: E8385JA (with gel) E8385CE (without gel) GE Medical Systems 800-558-5102
Civco Sterile Ultrasound Transducer Cover Kit	<ul style="list-style-type: none"> • 5.5" width x 58" length • Contains sterile gel packet • Latex-free 	Part #: 610-575 Civco Medical Instruments Kalona, Iowa 800-445-6741
Ultrasound Gel		
Sterile Gel Packets	<ul style="list-style-type: none"> • 20 ml packets, 48/box 	Part #: 610-117 Civco Medical Instruments Kalona, Iowa 800-445-6741
Non-sterile gel	<ul style="list-style-type: none"> • .25 liter bottles, 12/case • 5 liters w/ 1.25 liter dispenser bottle 	Part #: 610-130 (.25L bottles) 610-116 (5L dispenser bottle) Civco Kalona, IA 800-445-6741
Pressure Bag		
IV Pressure Infusor Bag	<ul style="list-style-type: none"> • 1000ml • With pressure gauge 	Part #: 7789 Alliance Medical Russellville, MO 888-633-6908