

VarExTM

> Phlebectomy Instruments



Utilize the wide variety of VarEx instrument shapes and sizes to enhance current technique



Double-ended instruments:

These devices offer two tools in one. Flat spatula tip is especially helpful in dissecting tissue. (Similar to Varady style)



Blunt-tipped instruments:

These feature a straight shaft and are ideal for placing under a vein and lifting upwards. Blunt tip allows user to avoid unintentionally hooking dermis or other tissues. (Similar to Mueller style)



Sharp-tipped instruments:

These small, fine tips enable the user to pierce and exteriorize the vein, reducing trauma to surrounding tissues. Hook angle facilitates dissection. Ideal for smaller veins. (Similar to Ramelet style)



Crochet-style instruments:

The crochet instruments offer an alternate tip style to our smooth and sharp-tipped families.

-
- > **Convenience** – VNUS is your single source for high quality products to meet your needs in treating venous disease
 - > **Expertise** – Knowledgeable personnel and focused products provide solutions for your practice



VNUSTM
MEDICAL TECHNOLOGIES

VarExTM

> Phlebectomy Instruments

Benefits:

- High quality 440A stainless surgical grade steel
- Unique octagonal handle design permits user to roll the instrument between the fingers, reducing wrist stress
- Flat, ridged surfaces for ease of gripping

Ordering Information

TIP DESIGN	MODEL NUMBER
VarEx Crochet Hooks	
Blunt, Size 1 (smallest)	VPH-CH1
Blunt, Size 2 (largest)	VPH-CH2
VarEx SmoothTip Family	
Blunt, left-handed, Size 1 (smallest)	VPH-M1L
Blunt, left-handed, Size 2	VPH-M2L
Blunt, left-handed, Size 3	VPH-M3L
Blunt, left-handed, Size 4 (largest)	VPH-M4L
Blunt, right-handed, Size 1	VPH-M1R
Blunt, right-handed, Size 2	VPH-M2R
Blunt, right-handed, Size 3	VPH-M3R
Blunt, right-handed, Size 4	VPH-M4R
VarEx FineTip Family	
Sharp, left-handed, Size 1	VPH-R1L
Sharp, left-handed, Size 2	VPH-R2L
Sharp, right-handed, Size 1	VPH-R1R
Sharp, right-handed, Size 2	VPH-R2R
VarEx Duo Family	
Double-Ended, Dissector	VPH-VA1
Double-Ended, Extractor - small hook	VPH-VA2
Double-Ended, Extractor - large hook	VPH-VA3
Double-Ended, Microspatula	VPH-VA4
VarEx Forceps	
Mosquito Forceps – 5.5" curved	VPH-MOS1
Mosquito Forceps – 5.5" straight	VPH-MOS2



VarEx Crochet Tips



VarEx FineTips



VarEx Custom Grip



VarEx SmoothTips

> To order, call toll free: 1-888-797-VEIN (8346)



Corporate Office
 VNUS Medical Technologies
 5799 Fontanoso Way
 San Jose, CA 95138
 Tel: 408-360-7200
 Fax: 408-365-8480
 www.vnus.com